

COMPRESSION
THERAPY
COLLECTION



TRUFORM[®]
LEG HEALTH



CONTENTS

Compression Profile	3
Therapy Guide	4
Measuring	5

Fashion Wear

Ladies' Sheer.....	6
Ladies' Opaque	10
Ladies' Pattern Socks	14
Men's Dress Socks.....	16
Men's & Ladies' Casual.....	18
Unisex Travel Series	20

NEW! Microfiber Medical Leg Wear

Classic Medical	21
-----------------------	----

Medical Leg Wear

Classic Medical	24
Anti-Embolism.....	30

Specialty Wear

Diabetic & Health Socks	32
Upper Extremity.....	35

TRUFORM Accessories	36
---------------------------	----

TRUFORM Product Summary	38
-------------------------------	----

COMPRESSION PROFILE

TRUFORM'S 6 POINTS OF COMPRESSION

To be therapeutically effective, medical stockings must provide a pressure reduction from the ankle upwards.

Truform takes compression measurements at precise intervals between the ankle and proximal end of the stocking - 4 points for below knee, 6 points for thigh length and pantyhose styles - to insure an accurate distribution of the pressure gradient.

Multiple measurements guarantee that Truform garments are comfortable and will deliver the prescribed graduated elastic compression therapy, without undesirable pressure inversions, when properly sized and fit.



PROXIMAL THIGH ▶

MID THIGH ▶

FIBULAR HEAD ▶

MID CALF ▶

LOWER CALF ▶

ANKLE ▶



THERAPY GUIDE

There is an appropriate Truform elastic compression therapy product for each of your customers.

Mild: 8-15 mmHg

Leg Fatigue
Slight ankle, leg and foot swelling

Moderate: 15-20 mmHg

Minor varicosities
Minor varicosities during pregnancy
Tired legs
Minor ankle, leg and foot swelling
Post surgery

Anti-Embolism: 18 mmHg

Prevention of thromboembolism
Reclining or bed confined home convalescence

Firm: 20-30 mmHg

Tired, aching legs
Moderate varicosities
Moderate varicosities during pregnancy
Moderate edema
Post sclerotherapy and other surgery
Help prevent recurrence of venous ulcerations
Superficial thrombophlebitis

Extra Firm: 30-40 mmHg

Severe varicosities
Severe edema
Lymphatic edema
Helps prevent and manage venous ulcerations
Chronic venous insufficiency
Helps prevent post thrombotic syndrome
Manage manifestations of PTS
Orthostatic hypotension
Post sclerotherapy
Post surgery

COMPRESSION PACKAGING LABEL



AVAILABLE STYLES



Knee High



Thigh High




Pantyhose

All compression ranges stated are the mean compression for an average ankle for a specific size. Please refer to the product insert, product label and package for comprehensive instructions on the proper selection and safe use of Truform products. Sizing for Truform support hosiery is determined by combinations of circumference and length measurements, or by combinations of height and weight data. Some items, due to their design, may be "quick-fit" by shoe or sock size. Sizing charts for Truform hosiery are located on each product page throughout this catalog as appropriate.

MEASURING


Please follow these instructions carefully to determine the correct size for different product styles. Measurements should always be taken in the following order, and as soon as possible after waking when swelling is at its minimum.



1 Ankle circumference



2 Calf circumference



3 Knee length from the right-angle bend of the knee to the floor



4 Widest circumference of the thigh



5 Leg length from the floor up to the mid-thigh (1) or to the gluteal fold (2) for Anti-Embolism products.

Scan this QR Code to view the TRUFORM fitting video:

TRUFORM
LEG HEALTH

TRUFORM
LEG HEALTH

GENERAL GUIDELINES

- Measurements should always be taken as soon as possible after waking when swelling is at its minimum
- Measure against bare skin
- The measuring system is based on having the foot flat against the floor and leg at a right angle to the foot
- Measurements must be taken by pulling the tape snugly around the leg but without constricting the skin

LADIES' SHEER

1775 15-20 mmHg
Pantyhose BLACK



LADIES' SHEER KNEE HIGH

CLOSED TOE

- 1763 LITES 8-15 mmHg
- 1773 LITES 15-20 mmHg
- 1782 Dot Pattern 15-20 mmHg
- 1783 Diamond Pattern 15-20 mmHg
- 0263 TruSheer 20-30 mmHg
- 0253 TruSheer 30-40 mmHg

OPEN TOE

- 1762 LITES 8-15 mmHg
- 1772 LITES 15-20 mmHg



1782 DOT PATTERN



1783 DIAMOND PATTERN



1763 BLACK



1772 NUDE

Material Information

1763
88% nylon and 12% spandex

1772 | 1773 | 1782 | 1783
80% nylon and 20% spandex

0263
78% nylon and 22% spandex

0253
57% nylon and 43% spandex

Product Description

Provide mild, moderate, firm, or extra firm graduated compression to help energize the legs, and reduce leg discomfort. They provide long-lasting support that appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Mild, Moderate, Firm, or Extra-Firm support
- Fashionable sheer fabric
- Knit-in heel (1772, 1773, 1782, 1783, 0253, 0263)
- Available in diamond or dot patterns (1782, 1783)

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

LATEX FREE

	SHOE SIZE	
SMALL	4 - 5	1762 1763
MEDIUM	5 ½ - 7 ½	
LARGE	8 - 10 ½	
X-LARGE	11 - 12 ½	
2X-LARGE	13 and over	

LADIES' SHEER KNEE HIGH COMPRESSION STOCKINGS			
MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
1762	8 - 15 mmHg	Knee High Open Toe	✓
1763	8 - 15 mmHg	Knee High	✓
1772	15 - 20 mmHg	Knee High Open Toe	✓
1773	15 - 20 mmHg	Knee High	✓
1782	15 - 20 mmHg	Knee High Dot Pattern	✓
1783	15 - 20 mmHg	Knee High Diamond Patten	✓
0263	20 - 30 mmHg	Knee High	✓
0253	30 - 40 mmHg	Knee High	✓

	CIRCUMFERENCE		
	ANKLE	CALF	
SMALL	7" - 8 ¼"	11" - 15"	0263 0253 1772 1773 1782 1783
MEDIUM	8 ¼" - 9 ½"	11 ¾" - 16 ½"	
LARGE	9 ½" - 10 ¾"	12 ½" - 18"	
X-LARGE	10 ¾" - 12"	13" - 19"	


	CIRCUMFERENCE		
	ANKLE	CALF	
2X-LARGE	12 ½" - 13 ¾"	19 ½" - 22"	1772 1773
3X-LARGE	12 ½" - 13 ¾"	21 ½" - 24"	

LADIES' SHEER THIGH HIGH

CLOSED TOE

- 1764 LITES 8-15 mmHg
- 1774 LITES 15-20 mmHg
- 0264 TruSheer 20-30 mmHg
- 0254 TruSheer 30-40 mmHg

Silicone stay-up top



264 BLACK

Material Information

1764
88% nylon and 12% spandex

1774
83% nylon and 17% spandex

0264 | 0254
57% nylon and 43% spandex

Product Description


Provide mild, moderate, firm, or extra firm graduated compression to help energize the legs, and reduce leg discomfort. They provide long-lasting support that appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Mild, Moderate, Firm, or Extra-Firm support
- Fashionable sheer fabric
- Attractive lace stay-up top

LATEX FREE

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

LADIES' SHEER THIGH HIGH COMPRESSION STOCKINGS			
MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
1764	8 - 15 mmHg	Thigh High	✓
1774	15 - 20 mmHg	Thigh High	✓
0264	20 - 30 mmHg	Thigh High	✓
0254	30 - 40 mmHg	Thigh High	✓

CIRCUMFERENCE

	ANKLE	CALF	THIGH	
SMALL	7" - 8 1/4"	11" - 15"	15 3/4" - 24 1/2"	0264 0254
MEDIUM	8 1/4" - 9 1/2"	11 3/4" - 16 1/2"	18" - 27 1/2"	
LARGE	9 1/2" - 10 3/4"	12 1/2" - 18"	21 1/4" - 30 3/4"	
X-LARGE	10 3/4" - 12"	13" - 19"	23 1/2" - 32"	

CIRCUMFERENCE

	ANKLE	CALF	THIGH	
SMALL	7" - 8 1/4"	11" - 14"	up to 18 1/2"	1774
MEDIUM	8 1/4" - 9 3/4"	13 1/2" - 16"	up to 19 1/2"	
LARGE	9 3/4" - 11"	15 1/2" - 18"	up to 20 1/2"	
X-LARGE	11" - 12 1/2"	17 1/2" - 20"	up to 21 1/2"	

CIRCUMFERENCE

	SHOE SIZE	THIGH	
SMALL	4 - 5	18 1/2"	1764
MEDIUM	5 1/2 - 7 1/2	19 1/2"	
LARGE	8 - 10 1/2	20 1/2"	
X-LARGE	11 - 12 1/2	21 1/2"	

LADIES' SHEER PANTYHOSE

CLOSED TOE

- 1765 LITES 8-15 mmHg
- 1775 LITES 15-20 mmHg
- 0265 TruSheer 20-30 mmHg
- 0255 TruSheer 30-40 mmHg

MATERNITY FIGURE

- 1777 LITES 15-20 mmHg
- 0267 TruSheer 20-30 mmHg



1765 BLACK

Material Information

1765
88% nylon and 12% spandex

1775 | 1777
84% nylon and 16% spandex

0267
71% nylon and 29% spandex

0265 | 0255
57% nylon and 43% spandex


Product Description

Provide mild, moderate, firm, or extra firm graduated compression to help energize the legs, and reduce leg discomfort. They provide long-lasting support that appeals to today's active, fashion conscious woman - excellent for business or casual wear.


Product Features

- Therapeutic graduated compression
- Mild, Moderate, Firm, or Extra-Firm support
- Fashionable sheer fabric
- Flat panty seams for a comfortable inconspicuous fit
- Ventilated cotton panel for comfort
- Available with maternity panel (1777 and 0267)

AVAILABLE COLORS *



*Colors may not be available in all styles.
Please refer to the Product Summary on pages 38 and 39 for more information.



LADIES' SHEER COMPRESSION PANTYHOSE			
MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
1765	8 - 15 mmHg	Pantyhose	✓
1775	15 - 20 mmHg	Pantyhose	✓
1777	15 - 20 mmHg	Maternity Pantyhose	✓
0265	20 - 30 mmHg	Pantyhose	✓
0255	30 - 40 mmHg	Pantyhose	✓
0267	20 - 30 mmHg	Maternity Pantyhose	✓

	CIRCUMFERENCE		
	ANKLE	THIGH	HIPS
SMALL	7" - 8 1/4"	15 3/4" - 24 1/2"	28" - 46"
MEDIUM	8 1/4" - 9 1/2"	18" - 27 1/2"	30" - 50"
LARGE	9 1/2" - 10 3/4"	21 1/2" - 30 3/4"	32" - 54"
X-LARGE	10 3/4" - 12"	23 1/2" - 32"	40" - 65"

LADIES' OPAQUE

0373 15-20 mmHg
Knee High BLACK



LADIES' OPAQUE KNEE HIGH

CLOSED TOE

- 0373 15-20 mmHg
- 0363 20-30 mmHg

OPEN TOE

- 0371 15-20 mmHg
- 0361 20-30 mmHg



0373 BLACK

0361 BEIGE

Material Information

0371 | 0373 | 0361 | 0363
85% nylon and 15% spandex


Product Description

TRUFORM Opaque resemble high fashion hosiery yet feature moderate and firm graduated compression technology - very soft for the most comfortable stocking in its compression class. The "easy-on", "easy-off" construction is lightweight and highly effective in helping restore leg health. Long-lasting medical-grade support, fashionable look and super soft texture appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Moderate or Firm support
- Advanced microfiber-soft yarns and knit-in heel and toe offer comfort, better fit and durability
- Fashionably smooth and soft opaque fabric provides for year-round wear and covers blemishes
- Soft, stay-up top for all-day comfort

AVAILABLE COLORS



CIRCUMFERENCE

	ANKLE	CALF	
SMALL	7" - 8 1/4"	11" - 15"	0373 0363 0371 0361
MEDIUM	8 1/4" - 9 1/2"	11 3/4" - 16 1/2"	
LARGE	9 1/2" - 10 3/4"	12 1/2" - 18"	
X-LARGE	10 3/4" - 12"	13" - 19"	

LATEX FREE

LADIES' OPAQUE KNEE HIGH COMPRESSION STOCKINGS

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
0371	15 - 20 mmHg	Knee High Open Toe	✓
0373	15 - 20 mmHg	Knee High	✓
0361	20 - 30 mmHg	Knee High Open Toe	✓
0363	20 - 30 mmHg	Knee High	✓

LADIES' OPAQUE THIGH HIGH

CLOSED TOE

- 0374 15-20 mmHg
- 0364 20-30 mmHg

OPEN TOE

- 0362 20-30 mmHg

Silicone stay-up top



0364 BEIGE

0362 BLACK

Material Information

0364 | 0374 | 0362
85% nylon and 15% spandex

Product Description


TRUFORM Opaque resemble high fashion hosiery yet feature firm graduated compression technology - very soft for the most comfortable stocking in its compression class. The "easy-on", "easy-off" construction is lightweight and highly effective in helping restore leg health. Long-lasting medical-grade support, fashionable look and super soft texture appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Firm support
- Advanced microfiber-soft yarns and knit-in heel offer comfort, better fit and durability
- Fashionably smooth soft opaque fabric provides for year-round wear and helps conceal blemishes
- Functional silicone dot stay-up top for all-day wear

LATEX FREE

AVAILABLE COLORS



CIRCUMFERENCE

	ANKLE	CALF	THIGH	
SMALL	7" - 8 1/4"	11" - 15"	15 3/4" - 24 1/2"	0364 0374 0362
MEDIUM	8 1/4" - 9 1/2"	11 3/4" - 16 1/2"	18" - 27 1/2"	
LARGE	9 1/2" - 10 3/4"	12 1/2" - 18"	21 1/4" - 30 3/4"	
X-LARGE	10 3/4" - 12"	13" - 19"	23 1/2" - 32"	

LADIES' OPAQUE THIGH HIGH COMPRESSION STOCKINGS

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
0374	15 - 20 mmHg	Thigh High	✓
0364	20 - 30 mmHg	Thigh High	✓
0362	20 - 30 mmHg	Thigh High Open Toe	✓

LADIES' OPAQUE PANTYHOSE

CLOSED TOE

- 0375 15-20 mmHg
- 0365 20-30 mmHg



0365 BLACK

0365 BIEGE

Material Information

0365 | 0375
85% nylon and 15% spandex

Product Description


TRUFORM Opaque resemble high fashion hosiery yet feature firm graduated compression technology - very soft for the most comfortable stocking in its compression class. The "easy-on", "easy-off" construction is lightweight and highly effective in helping restore leg health. Long-lasting medical-grade support, fashionable look and super soft texture appeals to today's active, fashion conscious woman - excellent for business or casual wear.

Product Features

- Therapeutic graduated compression
- Firm support
- Advanced microfiber-soft yarns and knit-in heel and toe offer comfort, better fit and durability
- Fashionably smooth soft opaque fabric provides for year-round wear and helps conceal blemishes
- Flat panty seams for a comfortable inconspicuous fit
- Ventilated cotton panel for greater wearing comfort

LATEX FREE

AVAILABLE COLORS



CIRCUMFERENCE

	ANKLE	THIGH	HIPS
SMALL	7" - 8 1/4"	15 3/4" - 24 1/2"	28" - 46"
MEDIUM	8 1/4" - 9 1/2"	18" - 27 1/2"	30" - 50"
LARGE	9 1/2" - 10 3/4"	21 1/4" - 30 3/4"	32" - 54"
X-LARGE	10 3/4" - 12"	23 1/2" - 32"	40" - 65"

0365
0375

LADIES' OPAQUE COMPRESSION PANTYHOSE

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
0375	15 - 20 mmHg	Pantyhose	✓
0365	20 - 30 mmHg	Pantyhose	✓

LADIES' PATTERN SOCKS



1976 15-20 mmHg
Knee High Diamond
BLACK



LADIES' PATTERN SOCKS KNEE HIGH

RIB

● 1973 15-20 mmHg

CABLE

● 1975 15-20 mmHg

DIAMOND


● 1976 15-20 mmHg



1973 NAVY

1975 WHITE

1976 BLACK



Material Information

1973 | 1975 | 1976
85% nylon and 15% spandex

Product Description


TRUFORM trouser socks are designed to help energize the legs and reduce leg discomfort. Microfiber nylon on the rib, cable, and diamond patterns is soft, and cooler to wear than standard nylon. Attractive designer-knit delicate patterns make these socks an excellent choice for business, or business-casual wear.

Product Features

- Therapeutic graduated compression
- Moderate support
- Soft, comfortable microfiber nylon construction
- Smooth, seamless linked toe closure
- Extra-wide comfort band
- Attractive, designer-knit delicate patterns

LATEX FREE

AVAILABLE COLORS



CIRCUMFERENCE

	ANKLE	CALF	
SMALL	7" - 8 1/4"	11" - 14"	1973 1975 1976
MEDIUM	8 1/4" - 9 3/4"	13 1/2" - 16"	
LARGE	9 1/2" - 11"	15 1/2" - 18"	
X-LARGE	11" - 12 1/2"	17 1/2" - 20"	

LADIES' KNEE HIGH COMPRESSION SOCKS

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
1973	15 - 20 mmHg	Knee High Rib	✓
1975	15 - 20 mmHg	Knee High Cable	✓
1976	15 - 20 mmHg	Knee High Diamond	✓

MEN'S DRESS



1944 20-30 mmHg
Knee High BLACK

MEN'S DRESS KNEE | THIGH HIGH

KNEE HIGH

- 1942 8-15 mmHg
- 1943 15-20 mmHg
- 1944 20-30 mmHg
- 1954 30-40 mmHg

THIGH HIGH

- 1945 20-30 mmHg



Silicone stay-up top



Material Information

1942
93% nylon and 7% spandex

1943
90% nylon and 10% spandex

1944 | 1954 | 1945
86% nylon and 14% spandex

Product Description


Truform Dress-style Support Socks for men provide mild, moderate, firm or extra firm graduated compression to help energize the legs - ideal when sitting or standing for long periods. Nylon-spandex materials are soft, and resemble specialty store, designer-knit socks. Contemporary rib-knit styling is suitable for business or business-casual settings.

Product Features

- Therapeutic graduated compression
- Mild, Moderate, Firm or Extra firm support available
- Roomy toe and defined heel for comfort and correct fit
- Contemporary rib-knit styling
- Functional, silicone dot stay-up top for all-day wear (1945)
- **NEW! Model 1943** now made using enhanced microfiber material resulting in improved stretch and fit, lightweight feel, ease of application and wearing comfort.

LATEX FREE

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

MEN'S KNEE/THIGH HIGH COMPRESSION SOCKS

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
1942*	8 - 15 mmHg	Knee High	✓
1943	15 - 20 mmHg	Microfiber Knee High	✓
1944	20 - 30 mmHg	Knee High	✓
1945	20 - 30 mmHg	Thigh High	✓
1954	30 - 40 mmHg	Knee High	✓

CIRCUMFERENCE

	SHOE SIZE	ANKLE	CALF	MID-THIGH	LENGTH TO KNEE	
SMALL	7 and under	7" - 8 1/4"	11" - 14"	15 3/4" - 24 1/2"	up to 15"	1942
MEDIUM	7 1/2 - 10	8 1/4" - 9 3/4"	13 1/2" - 16"	18" - 27 1/2"	up to 16"	1943
LARGE	10 1/2 - 12	9 3/4" - 11"	15 1/2" - 18"	21 1/2" - 30 3/4"	up to 17"	1944
X-LARGE	12 and over	11" - 12 1/2"	17 1/2" - 20"	23 1/2" - 32"	up to 18"	1945 1954

*Use shoe size for 1942

MEN'S AND LADIES' CASUAL



1933 15-20 mmHg
Casual Style Knee High BLACK



MEN'S & LADIES' CASUAL SOCKS

LADIES' KNEE HIGH

● **1963** 10-20 mmHg

MEN'S KNEE HIGH

● **1933** 15-20 mmHg

White color option available in athletic box

● **1934** 20-30 mmHg

MEN'S CREW

● **1932** 15-20 mmHg



1963 NAVY

1933 TAN

1932 WHITE

Material Information

1932 | 1933
64% acrylic, 29% nylon and 7% spandex

1934
60% polyester, 31% nylon and 9% spandex

1963
60% acrylic, 33% nylon and 7% spandex

Product Description


Truform Men's & Ladies' Casual Support Socks are designed to help improve circulation and reduce leg discomfort. Shock-absorbing cushioned toe, sole and heel, moisture wicking fibers and anti-microbial, stay-fresh finish provide all-day wearing comfort. Excellent for court sports, golfing, jogging, walking, or any leisure activity that requires comfort and support.

Product Features

- Therapeutic graduated compression
- Moderate support
- Shock-absorbing cushioned toe, sole and heel
- Moisture wicking acrylic materials
- Anti-microbial, stay-fresh treatment helps control odor

LATEX FREE

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

MEN'S AND LADIES' CASUAL COMPRESSION SOCKS

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
1963 Ladies'	10 - 20 mmHg	Casual Style Knee High	✓
1932 Men's (Colors)	15 - 20 mmHg	Casual Style Crew Length	✓
1933 Men's (White)	15 - 20 mmHg	Athletic Style Crew Length	✓
1933 Men's (Colors)	15 - 20 mmHg	Casual Style Knee High	✓
1933 Men's (White)	15 - 20 mmHg	Athletic Style Knee High	✓
1934 Men's	20 - 30 mmHg	Casual Style Knee High	✓

CIRCUMFERENCE

	ANKLE	CALF
SMALL	7" - 8 ¼"	11" - 14"
MEDIUM	8 ¼" - 9 ¾"	13 ½" - 16"
LARGE	9 ¾" - 11"	15 ½" - 18"
X-LARGE	11" - 12 ½"	17 ½" - 20"

1932
1933
1934

CIRCUMFERENCE

	ANKLE	CALF
SMALL	7" - 8 ¼"	11" - 15"
MEDIUM	8 ¼" - 9 ½"	11 ¾" - 16 ½"
LARGE	9 ½" - 10 ¾"	12 ½" - 18"
X-LARGE	10 ¾" - 12"	13" - 19"

1963

TRAVEL SERIES



TRUFORM[®]
TRAVEL SERIES

UNISEX TRAVEL SERIES KNEE HIGH

KNEE HIGH

● 1923 15-20 mmHg



1923 BROWN



1923 CHARCOAL

Material Information

1923
90% nylon and 10% spandex

Product Description

TRUFORM® TRAVEL SERIES graduated compression socks are specifically designed for travel. They were developed using a design that helps improve circulation and reduces swelling while in close or cramped quarters for extended periods of time. Tested by pilots and people of all ages traveling worldwide, TRUFORM® is sure to be your ideal travel companion!

Product Features

- Therapeutic graduated compression
- Moderate support
- Roomy toe and defined heel for comfort and correct fit
- Contemporary rib-knit styling

LATEX FREE

AVAILABLE COLORS

Charcoal
CH



Brown
BN



Navy
NV



Black
BL



UNISEX KNEE HIGH TRAVEL COMPRESSION SOCKS

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
1923	15 - 20 mmHg	Knee High	✓

	CIRCUMFERENCE		SHOE SIZE			
	ANKLE	CALF	US MEN'S	US WOMEN'S	UK	EURO
SMALL	7" - 8 ¼"	11" - 14"	4 ½ - 7	6 - 8 ½	3 ½ - 6	36 - 40
MEDIUM	8 ¼" - 9 ¾"	13 ½" - 16"	7 ½ - 10 ½	9 - 11 ½	6 ½ - 9	40 ½ - 43
LARGE	9 ¾" - 11"	15 ½" - 18"	10 ½ - 12	12 - 13 ½	9½ - 11	43 ½ - 45
X-LARGE	11" - 12 ½"	17 ½" - 20"	12 ½ & over	14 & over	11 ½ & over	45 ½ & over

MICROFIBER MEDICAL LEG WEAR



NEW FROM TRUFORM! **MICROFIBER MEDICAL** Knee High 20-30 mmHg Compression Socks



**ENHANCED
STRETCH & FIT**



**SUPERIOR
LIGHTWEIGHT FEEL**



**IMPROVED
EASE OF APPLICATION**



**ULTIMATE
WEARING COMFORT**



MICROFIBER MEDICAL LEG WEAR BELOW KNEE

BELOW KNEE CLOSED TOE SOFT TOP

● 9965 20-30 mmHg

BELOW KNEE OPEN TOE SOFT TOP

● 0965 20-30 mmHg



0965 BEIGE

9965 BLACK

Material Information

9965 | 0965
84% nylon and 16% spandex

Product Description

Truform MicroFiber Medical Compression Socks utilize graduated compression technology to help restore leg health – ideal when sitting or standing for long periods.


These socks are made using premium four-way stretch air-covered MicroFiber yarns, which combine elastic with non-elastic materials that significantly improve appearance, function and durability.

Product Features

- Therapeutic graduated compression
- Firm support
- Enhanced stretch and fit
- Superior lightweight feel
- Improved ease of application
- Air-covered Microfiber yarns provide the ultimate wearing comfort

LATEX FREE

AVAILABLE COLORS



CIRCUMFERENCE

	ANKLE	CALF
SMALL	7" - 8 ¼"	11" - 14"
MEDIUM	8 ¾" - 9 ⅝"	13 ½" - 16"
LARGE	9¾" - 11"	15 ½" - 18"
X-LARGE	11⅞" - 12 ¾"	17 ½" - 20"
2X-LARGE	12 ½" - 13 ¾"	19 ½" - 22"
3X-LARGE	12 ½" - 13 ¾"	21 ½" - 24"

MICROFIBER MEDICAL KNEE HIGH COMPRESSION SOCKS

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
9965	20 - 30 mmHg	Below Knee Closed Toe	✓
0965	20 - 30 mmHg	Below Knee Open Toe	✓

MEDICAL LEG WEAR



864 20-30 mmHg
Knee High Open Toe
BEIGE



MEDICAL LEG WEAR BELOW KNEE CLOSED TOE

BELOW KNEE SOFT TOP

- 8875 15-20 mmHg
- 8865 20-30 mmHg
- 8865S 20-30 mmHg
- 8845 30-40 mmHg
- 8845S 30-40 mmHg

BELOW KNEE DOT TOP

- 8864 20-30 mmHg
- 8844 30-40 mmHg

Silicone dot top



8865 BEIGE



Material Information

8875
85% nylon and 15% spandex

8865 | 8865S | 8845 | 8845S | 8864 | 8844
80% nylon and 20% spandex

Product Description

Designed to help relieve moderate and pronounced conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm or Extra Firm support
- Classic, concealing smooth knit texture fabric
- Available in soft top and non-constricting, silicone dot top
- Standard and larger sizes

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

LATEX FREE

	CIRCUMFERENCE		LENGTH	
	ANKLE	CALF	TO KNEE	
SMALL	7" - 8 1/4"	11" - 14"	<p>Regular: Greater than 15" Short: 15" or less</p>	
MEDIUM	8 3/8" - 9 5/8"	13 1/2" - 16"		
LARGE	9 3/4" - 11"	15 1/2" - 18"		
X-LARGE	11 1/8" - 12 3/8"	17 1/2" - 20"		
2X-LARGE	12 1/2" - 13 3/4"	19 1/2" - 22"		up to 18"
3X-LARGE	12 1/2" - 13 3/4"	21 1/2" - 24"		up to 18"
4X-LARGE	13" - 14 1/4"	23 1/2" - 26"	up to 18"	

*8875 not available in 2XL, 3XL and 4XL

CLASSIC MEDICAL KNEE HIGH COMPRESSION SOCKS			
MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
8875	15 - 20 mmHg	Below Knee Soft Top	✓
8864	20 - 30 mmHg	Below Knee Dot Top	✓
8865	20 - 30 mmHg	Below Knee Soft Top	✓
8865S (Short)	20 - 30 mmHg	Below Knee Soft Top	✓
8844	30 - 40 mmHg	Below Knee Dot Top	✓
8845	30 - 40 mmHg	Below Knee Soft Top	✓
8845S (Short)	30 - 40 mmHg	Below Knee Soft Top	✓

MEDICAL LEG WEAR BELOW KNEE OPEN TOE

BELOW KNEE OPEN TOE SOFT TOP

- **0875** 15-20 mmHg
- **0865** 20-30 mmHg
- **0865S** 20-30 mmHg
- **0845** 30-40 mmHg
- **0845S** 30-40 mmHg

BELOW KNEE OPEN TOE DOT TOP

- **0864** 20-30 mmHg
- **0844** 30-40 mmHg

Silicone dot top



0865 BEIGE

Material Information

0865 | 0865S | 0845 | 0845S | 0864 | 0844
80% nylon and 20% spandex

0875
85% nylon and 15% spandex

Product Description


Designed to help relieve moderate and pronounced conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. The open toe feature allows for easy toe inspection and permits one standard size to fit a variety of foot dimensions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm or Extra Firm support
- Classic, concealing smooth knit texture fabric
- Available in soft top and non-constricting, silicone dot top
- Standard and larger sizes

LATEX FREE

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

CIRCUMFERENCE

LENGTH

	CIRCUMFERENCE		LENGTH
	ANKLE	CALF	TO KNEE
SMALL	7" - 8 1/4"	11" - 14"	 Regular: Greater than 15" Short: 15" or less
MEDIUM	8 3/8" - 9 5/8"	13 1/2" - 16"	
LARGE	9 3/4" - 11"	15 1/2" - 18"	up to 18"
X-LARGE	11 1/8" - 12 3/8"	17 1/2" - 20"	
2X-LARGE	12 1/2" - 13 3/4"	19 1/2" - 22"	up to 18"
3X-LARGE	12 1/2" - 13 3/4"	21 1/2" - 24"	
4X-LARGE	13" - 14 1/4"	23 1/2" - 26"	up to 18"

*0875 not available in 2XL, 3XL and 4XL

CLASSIC MEDICAL KNEE HIGH OPEN TOE COMPRESSION SOCKS

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
0875	15 - 20 mmHg	Below Knee Soft Top	✓
0864	20 - 30 mmHg	Below Knee Dot Top	✓
0865	20 - 30 mmHg	Below Knee Soft Top	✓
0865S (Short)	20 - 30 mmHg	Below Knee Soft Top	✓
0844	30 - 40 mmHg	Below Knee Dot Top	✓
0845	30 - 40 mmHg	Below Knee Soft Top	✓
0845S (Short)	30 - 40 mmHg	Below Knee Soft Top	✓

MEDICAL LEG WEAR THIGH HIGH CLOSED TOE

THIGH HIGH SOFT TOP

- 8866 20-30 mmHg
- 8846 30-40 mmHg

THIGH HIGH DOT TOP

- 8868 20-30 mmHg
- 8848 30-40 mmHg

THIGH HIGH LACE TOP

- 8867 20-30 mmHg

Silicone lace
Stay-up top



Silicone dot
top



8867 BEIGE

8868 BEIGE



Material Information

8866 | 8868 | 8867
80% nylon and 20% spandex

8846 | 8848
83% nylon and 17% spandex

Product Description


Designed to help relieve moderate and pronounced conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm or extra firm support
- Classic, concealing smooth knit texture fabric
- Available in soft top and non-constricting, silicone dot top

Note: A garter belt should be worn with soft top style numbers 8846 and 8866

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

LATEX FREE

	CIRCUMFERENCE			LENGTH
	ANKLE	CALF	THIGH	TO THIGH
SMALL	7" - 8 ¼"	11" - 14"	up to 20"	up to 26"
MEDIUM	8 ⅜" - 9 ⅝"	13 ½" - 16"	up to 22"	up to 28"
LARGE	9 ¾" - 11"	15 ½" - 18"	up to 23"	up to 29"
X-LARGE	11 ⅛" - 12 ⅜"	17 ½" - 20"	up to 24"	up to 31"
2X-LARGE	12 ½" - 13 ½"	21 ½" - 24"	up to 32"	up to 33"

CLASSIC MEDICAL THIGH HIGH COMPRESSION SOCKS			
MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
8866	20 - 30 mmHg	Thigh High Soft Top	✓
8867	20 - 30 mmHg	Thigh High Lace Top	✓
8868	20 - 30 mmHg	Thigh High Dot Top	✓
8846	30 - 40 mmHg	Thigh High Soft Top	✓
8848	30 - 40 mmHg	Thigh High Dot Top	✓

MEDICAL LEG WEAR THIGH HIGH OPEN TOE

THIGH HIGH SOFT TOP

- 0866 20-30 mmHg
- 0846 30-40 mmHg

THIGH HIGH DOT TOP

- 0868 20-30 mmHg
- 0848 30-40 mmHg

Silicone
dot top



0868 BEIGE



Material Information

0866 | 0868 | 0846 | 0848
83% nylon and 17% spandex

Product Description

Designed to help relieve moderate and pronounced conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. The open toe feature allows for easy toe inspection and permits one standard size to fit a variety of foot dimensions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm or extra firm support
- Classic, concealing smooth knit texture fabric
- Available in soft top and non-constricting, silicone dot top

Note: A garter belt should be worn with soft top style numbers 0846 and 0866

LATEX FREE

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

	CIRCUMFERENCE			LENGTH
	ANKLE	CALF	THIGH	TO THIGH
SMALL	7" - 8 ¼"	11" - 14"	up to 20"	up to 26"
MEDIUM	8 ¾" - 9 ⅝"	13 ½" - 16"	up to 22"	up to 28"
LARGE	9 ¾" - 11"	15 ½" - 18"	up to 23"	up to 29"
X-LARGE	11 ⅛" - 12 ⅜"	17 ½" - 20"	up to 24"	up to 31"
2X-LARGE	12 ½" - 13 ½"	21 ½" - 24"	up to 32"	up to 33"

CLASSIC MEDICAL THIGH HIGH OPEN TOE COMPRESSION SOCKS

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
0866	20 - 30 mmHg	Thigh High Soft Top	✓
0868	20 - 30 mmHg	Thigh High Dot Top	✓
0846	30 - 40 mmHg	Thigh High Soft Top	✓
0848	30 - 40 mmHg	Thigh High Dot Top	✓

MEDICAL LEG WEAR PANTYHOSE

STANDARD FIGURE

● 1756 20-30 mmHg

MATERNITY FIGURE

● 1757 20-30 mmHg



FULL FIGURE

● 1758 20-30 mmHg



1756 BLACK

Material Information

1756 | 1757 | 1758
80% nylon and 20% spandex

Product Description

Designed to help relieve moderate conditions associated with poor circulation. Truform medical stockings feature classic styling and a smooth knit fabric that helps conceal unsightly conditions. Color options are basic neutrals for all-occasion wear. For men and women.

Product Features

- Therapeutic graduated compression
- Firm support.
- Classic, concealing smooth-knit texture fabric
- Flat panty seams for a comfortable inconspicuous fit
- Ventilated cotton panel for comfort
- Available with maternity panel (1757)

LATEX FREE



AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

CLASSIC MEDICAL COMPRESSION PANTYHOSE			
MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
1756	20 - 30 mmHg	Standard Figure	✓
1757	20 - 30 mmHg	Maternity Figure	✓
1758	20 - 30 mmHg	Full Figure	✓



ANTI-EMBOLISM




ANTI-EMBOLISM BELOW KNEE | THIGH LENGTH

BELOW KNEE

- 8808 | 8808S (Short) 18 mmHg

BELOW KNEE OPEN TOE

- 0808 | 0808S (Short) 18 mmHg



THIGH LENGTH

- 8810 18 mmHg

THIGH LENGTH OPEN TOE

- 0810 18 mmHg



Material Information

8808 | 8808S | 0808 | 0808S
89% nylon and 11% spandex

8810
85% nylon and 15% spandex

0810
82% nylon and 18% spandex

Product Description

Doctor-recommended for the recuperating patient, Truform anti-embolism Stockings are designed to help reduce the danger of thromboembolic disease for the reclining or bed-confined patient. They are also ideal for use with the post-surgical, home convalescing patient. The open toe feature allows for easy toe inspection by caregiver. Beige and black color options allow for a less clinical appearance during early ambulation.

Product Features

- Therapeutic graduated compression; 18 mmHg targeted at the ankle
- Non-constricting, soft stay-up top for comfort (8808, 0808)
- Open toe or closed toe designs
- Regular and Short lengths (Greater or less than 15" from heel to bend of knee)
- Pressure-relief panel in the upper thigh (8810, 0810)
- Silicone stay-up top (8810, 0810)

	CIRCUMFERENCE			LENGTH	
	ANKLE	CALF	THIGH	TO KNEE	TO THIGH
SMALL	7" - 8 ¼"	11" - 14"	up to 21"	 Regular: Greater than 15" Short: 15" or less	up to 28"
MEDIUM	8 ¾" - 9 ⅝"	13 ½" - 16"	up to 22"		up to 29"
LARGE	9 ¾" - 11"	15 ½" - 18"	up to 24"		up to 30"
X-LARGE	11 ⅙" - 12 ⅜"	17 ½" - 20"	up to 26"		up to 32"
2X-LARGE	12 ½" - 13 ¾"	19 ½" - 22"	N/A	up to 18"	N/A
3X-LARGE	12 ½" - 13 ¾"	21 ½" - 24"	N/A	up to 18"	N/A
4X-LARGE	13" - 14 ¼"	23 ½" - 26"	N/A	up to 18"	N/A

LATEX FREE

*0810 and 8810 not available in 2XL, 3XL and 4XL

* 0808 not available in 4XL

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

ANTI-EMBOLISM KNEE AND THIGH HIGH COMPRESSION			
MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
8808	18 mmHg	Below Knee	✓
8808S (Short)	18 mmHg	Below Knee Short	✓
0808	18 mmHg	Below Knee Open Toe	✓
0808S (Short)	18 mmHg	Below Knee Open Toe Short	✓
8810	18 mmHg	Thigh Length	✓
0810	18 mmHg	Thigh Length Open Toe	✓

DIABETIC AND HEALTH SOCKS



16-20-30 mmHg
Knee High Open Toe
Grey



TRUSOFT SOCKS



MINI-CREW 1911

CREW

● 1912 8-15 mmHg

KNEE HIGH

● 1913 8-15 mmHg

1913 BLACK

1912 BLACK

1911 WHITE

Material Information

1911 | 1912
81% acrylic 17% nylon and 2% spandex

1913
75% acrylic 22% nylon and 3% spandex

Product Description


TruSoft Socks surround the feet with full-cushioned comfort. Non-constricting, gentle, graduated compression helps keep sock in place to prevent bunching and sliding. Five gender-neutral sizes, moisture-wicking yarns, stay-fresh treatment, and NEW seamless linked toe closure help protect and comfort sensitive feet - ideal for the diabetic patient, or anyone with sensitive feet issues.

Product Features

- Ankle length, crew length, and knee high
- Mild graduated compression
- Moisture-wicking yarns
- Anti-microbial finish helps control odor
- NEW non-irritating seamless linked toe closure
- Full-cushioned foot

LATEX FREE

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

	MALE SHOE SIZE	FEMALE SHOE SIZE
X-SMALL	3 ½ - 5	4 ½ - 6 ½
SMALL	6 - 8	7 - 9
MEDIUM	8 ½ - 10	9 ½ - 11
LARGE	10 ½ - 12	11 ½ - 13
X-LARGE	12 ½ - 14	

1911
1912
1913

TRUSOFT DIABETIC SOCKS			
MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
1911	N/A	Mini Crew	✓
1912	8 - 15 mmHg	Crew Length	✓
1913	8 - 15 mmHg	Knee High	✓

COMFORT CARE DIABETIC SOCKS

LOOSE FIT TOP CREW

1918



1918 BLACK

Material Information

1918WH | 1918TN
80% cotton 15% polyester and 5% spandex

1918BL | 1918NV
80% cotton 15% nylon and 5% spandex

Product Description


Comfort Care casual socks are designed for anyone with special needs in foot care. Non-constricting, gentle, elasticity helps keep sock in place to prevent bunching and sliding. The natural cotton materials helps reduce irritation and blistering and provides exceptional comfort, easy care and excellent absorbency. Suitable for men and women - ideal for the diabetic patient, or anyone with sensitive feet issues.

Product Features

- Crew length, cut extra-wide and extra-deep
- Seamless toe closure
- Non-binding and comfortable
- Can be worn over a cast
- Suitable for wear at night for warmth
- Excellent cover for use with medicinal or softening creams

LATEX FREE

AVAILABLE COLORS



	MEN'S SHOE SIZE	LADIES' SHOE SIZE
SMALL	6.5 and under	5.5 and under
MEDIUM	6 - 9	7 - 10
LARGE	9 - 12	10 - 13
X-LARGE	12 ½ and over	13 ½ and over

COMFORT CARE DIABETIC SOCKS

MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
1918	N/A	Loose Fit Top Crew Length	✓

UPPER EXTREMITY

READY-WEAR SOFT TOP SLEEVE

- 3315 15-20 mmHg
- 3325 20-30 mmHg

READY-WEAR DOT TOP SLEEVE

- 3316 15-20 mmHg
- 3326 20-30 mmHg

READY-WEAR GAUNTLET

- 3312 15-20 mmHg
- 3322 20-30 mmHg



Material Information

3315 | 3325 | 3316 | 3326 | 3312 | 3322
73% nylon and 27% spandex

Product Description

Truform Ready To Wear Gauntlets and Sleeves are designed to help manage swelling following trauma, surgery, or radiation therapy. Lightweight microfiber knit construction allows for easy application, the soft air-permeable fabric provides extended wearing comfort. Available in three standard sizes and basic neutral colors.

Product Features

- Therapeutic, moderate and firm graduated compression
- Breathable, lightweight materials and construction is soft and easy to apply
- Silicone dot top to help sleeve remain in place
- For immediate wear with only 3 easy measurements
- Not made with natural rubber latex

MEASURING:

For gauntlet

Measure around the palm and around the smallest part of the wrist as shown.


For arm sleeve

Measure around the smallest part of the wrist, around the mid-lower arm (largest part of the forearm), and around the mid-upper arm as shown.



LATEX FREE

AVAILABLE COLORS *



*Colors may not be available in all styles.

Please refer to the Product Summary on pages 38 and 39 for more information.

READY-WEAR COMPRESSION UPPER EXTREMITY			
MODEL	COMPRESSION	DESCRIPTION	LATEX FREE
3315	15 - 20 mmHg	Arm Sleeve Soft Top	✓
3316	15 - 20 mmHg	Arm Sleeve Dot Top	✓
3325	20 - 30 mmHg	Arm Sleeve Soft Top	✓
3326	20 - 30 mmHg	Arm Sleeve Dot Top	✓
3312	15 - 20 mmHg	Gauntlet	✓
3322	20 - 30 mmHg	Gauntlet	✓

CIRCUMFERENCE

	PALM	WRIST	MID-LOWER	MID-UPPER	
SMALL	5 3/4" - 7 1/4"	5 1/4" - 7 1/2"	7" - 9"	8 1/2" - 15 1/4"	3315 3326 3325 3312 3316 3322
MEDIUM	7 1/2" - 8 1/2"	6" - 8 1/2"	9 1/4" - 11 1/4"	10 1/2" - 17 1/4"	
LARGE	8 3/4" - 10 1/4"	6 3/4" - 9 1/4"	11 1/4" - 13 1/4"	12 3/4" - 19 3/4"	



757 Donning Gloves

	PALM CIRCUMFERENCE
SMALL	6.5" - 7"
MEDIUM	7.5" - 8"
LARGE	8.5" - 9"
X-LARGE	9.5" - 10"

757



752 Stocking Applicator



750 Slip Sock

	REGULAR	TALL	LARGE
DONNING ZONE HEIGHT	9"	9"	10"
DONNING ZONE WIDTH	4"	4"	6"
ARM HEIGHT	13"	18"	18"
ARM DEPTH	3.5"	4"	4"

752

TRUFORM ACCESSORIES

750 SLIP SOCK

Often the most difficult part of applying medical compression stockings is getting them over the heel and ankle. The Slip Sock is designed to help open toe stockings slide smoothly onto the foot and over the ankle by reducing friction between the patient's skin and the compression stocking.

752 STOCKING APPLICATOR

The Truform stocking applicator allows donning of stockings with minimal bending and strain.

757 DONNING GLOVES

The non-slip, easy grip surface on the palm and fingers help ease compression hosiery on and protect the fabric from being damaged by rings or fingernails.

758 IT STAYS! ROLL-ON BODY ADHESIVE

A removable liquid adhesive with a handy roll-on applicator that is designed to hold support stockings and other garments in place. Non-irritating to skin, It Stays is ideal for medical stockings and hosiery, as well as other support products.



758 It Stays!

	WAIST CIRCUMFERENCE
X-SMALL	28" - 31"
SMALL	31" - 34"
MEDIUM	34" - 38"
MEDIUM/LARGE	38" - 42"
LARGE	42" - 46"
X-LARGE	46" - 50"

760
765

	WAIST CIRCUMFERENCE
SMALL	26" - 32"
MEDIUM	32" - 36"
LARGE	36" - 40"
X-LARGE	40" - 44"

762

TRUFORM GARTER BELTS

760 3" GARTER BELT

762 2" LIGHTWEIGHT GARTER BELT FOR WOMEN

765 6" GARTER BELT FOR MEN

These standard style garter belts are practical and easy to wear. Garter attachments are completely adjustable in length and provide a convenient way to keep thigh high stockings in place. Suitable for both men and women.



760 3" Garter Belt




765 6" Garter Belt for Men



762 2" Garter Belt



HIGH PRECISION MANUFACTURING



TRUFORM[®]

LEG HEALTH

KNITTING & SEWING OPERATIONS

Truform employs the latest methods of semi-automated manufacturing, and utilizes the best in technology to produce a world-class product at affordable prices - our prescribing physicians trust them, our many clients appreciate their extraordinary value.

Each machine is programmed to precise manufacturing specifications, and our technicians monitor every step of the knitting process. The combination of technology and knitting expertise ensures that Truform Stockings are among the highest quality, pressure gradient garments produced anywhere in the world.

We combine the finest materials and the most modern methods of 21st Century manufacturing, along with hand craftsmanship to produce Truform garments.

QUALITY ASSURANCE

Quality assurance checks are performed throughout our entire manufacturing process to verify fabric integrity, fit and a therapeutic, pressure-gradient profile.



TRUFORM PRODUCT SUMMARY

TRUFORM LITES SHEER SUPPORT HOSIERY PG. 6 - 9

MODEL	LENGTH	COMPRESSION	COLOR	SIZES
1762BG	Knee High	8 - 15 mmHg	Beige	Open Toe S-M-L-XL-2XL
1762BL	Knee High	8 - 15 mmHg	Black	Open Toe S-M-L-XL-2XL
1762NV	Knee High	8 - 15 mmHg	Navy	Open Toe S-M-L-XL-2XL
1762TP	Knee High	8 - 15 mmHg	Taupe	Open Toe S-M-L-XL-2XL
1763BG	Knee High	8 - 15 mmHg	Beige	S-M-L-XL-2XL
1763BL	Knee High	8 - 15 mmHg	Black	S-M-L-XL-2XL
1763NV	Knee High	8 - 15 mmHg	Navy	S-M-L-XL-2XL
1763TP	Knee High	8 - 15 mmHg	Taupe	S-M-L-XL-2XL
1764BG	Thigh High	8 - 15 mmHg	Beige	S-M-L-XL
1764BL	Thigh High	8 - 15 mmHg	Black	S-M-L-XL
1765BG	Pantyhose	8 - 15 mmHg	Beige	P-M-T-X-T-Q-QP
1765BL	Pantyhose	8 - 15 mmHg	Black	P-M-T-X-T-Q-QP
1765TP	Pantyhose	8 - 15 mmHg	Taupe	P-M-T-X-T-Q-QP
1772BG	Knee High	15 - 20 mmHg	Beige	Open Toe S-M-L-XL-2XL-3XL
1772IV	Knee High	15 - 20 mmHg	Ivory	Open Toe S-M-L-XL-2XL-3XL
1772LB	Knee High	15 - 20 mmHg	Light Beige	Open Toe S-M-L-XL-2XL-3XL
1772NV	Knee High	15 - 20 mmHg	Navy	Open Toe S-M-L-XL-2XL-3XL
1772TP	Knee High	15 - 20 mmHg	Taupe	Open Toe S-M-L-XL-2XL-3XL
1772ND	Knee High	15 - 20 mmHg	Nude	Open Toe S-M-L-XL-2XL-3XL
1772BL	Knee High	15 - 20 mmHg	Black	Open Toe S-M-L-XL-2XL-3XL
1772PK	Knee High	15 - 20 mmHg	Pink	Open Toe S-M-L-XL-2XL-3XL
1772TL	Knee High	15 - 20 mmHg	Teal	Open Toe S-M-L-XL-2XL-3XL
1772VT	Knee High	15 - 20 mmHg	Violet	Open Toe S-M-L-XL-2XL-3XL
1773BG	Knee High	15 - 20 mmHg	Beige	S-M-L-XL-2XL-3XL
1773BL	Knee High	15 - 20 mmHg	Black	S-M-L-XL-2XL-3XL
1773CH	Knee High	15 - 20 mmHg	Charcoal	S-M-L-XL-2XL-3XL
1773IV	Knee High	15 - 20 mmHg	Ivory	S-M-L-XL-2XL-3XL
1773LB	Knee High	15 - 20 mmHg	Light Beige	S-M-L-XL-2XL-3XL
1773ND	Knee High	15 - 20 mmHg	Nude	S-M-L-XL-2XL-3XL
1773NV	Knee High	15 - 20 mmHg	Navy	S-M-L-XL-2XL-3XL
1773TP	Knee High	15 - 20 mmHg	Taupe	S-M-L-XL-2XL-3XL
1773WH	Knee High	15 - 20 mmHg	White	S-M-L-XL-2XL-3XL
1773PK	Knee High	15 - 20 mmHg	Pink	S-M-L-XL-2XL-3XL
1773TL	Knee High	15 - 20 mmHg	Teal	S-M-L-XL-2XL-3XL
1773VT	Knee High	15 - 20 mmHg	Violet	S-M-L-XL-2XL-3XL
1774BG	Thigh High	15 - 20 mmHg	Beige	S-M-L-XL
1774BL	Thigh High	15 - 20 mmHg	Black	S-M-L-XL
1774CH	Thigh High	15 - 20 mmHg	Charcoal	S-M-L-XL
1774IV	Thigh High	15 - 20 mmHg	Ivory	S-M-L-XL
1774ND	Thigh High	15 - 20 mmHg	Nude	S-M-L-XL
1774TP	Thigh High	15 - 20 mmHg	Taupe	S-M-L-XL
1774WH	Thigh High	15 - 20 mmHg	White	S-M-L-XL
1775BG	Pantyhose	15 - 20 mmHg	Beige	P-M-T-X-T-Q-QP
1775BL	Pantyhose	15 - 20 mmHg	Black	P-M-T-X-T-Q-QP
1775IV	Pantyhose	15 - 20 mmHg	Ivory	P-M-T-X-T-Q-QP
1775LB	Pantyhose	15 - 20 mmHg	Light Beige	P-M-T-X-T-Q-QP
1775CH	Pantyhose	15 - 20 mmHg	Charcoal	P-M-T-X-T-Q-QP
1775ND	Pantyhose	15 - 20 mmHg	Nude	P-M-T-X-T-Q-QP
1775NV	Pantyhose	15 - 20 mmHg	Navy	P-M-T-X-T-Q-QP
1775TP	Pantyhose	15 - 20 mmHg	Taupe	P-M-T-X-T-Q-QP
1775WH	Pantyhose	15 - 20 mmHg	White	P-M-T-X-T-Q-QP
1777BG	Pantyhose	15 - 20 mmHg	Beige	Maternity P-M-T-X-T
1777BL	Pantyhose	15 - 20 mmHg	Black	Maternity P-M-T-X-T
1782BL	Knee High	15 - 20 mmHg	Black	Dot Pattern S-M-L-XL
1782CH	Knee High	15 - 20 mmHg	Charcoal	Dot Pattern S-M-L-XL
1782ES	Knee High	15 - 20 mmHg	Espresso	Dot Pattern S-M-L-XL
1782ND	Knee High	15 - 20 mmHg	Nude	Dot Pattern S-M-L-XL
1783BL	Knee High	15 - 20 mmHg	Black	Diamond Pattern S-M-L-XL
1783CH	Knee High	15 - 20 mmHg	Charcoal	Diamond Pattern S-M-L-XL
1783ES	Knee High	15 - 20 mmHg	Espresso	Diamond Pattern S-M-L-XL
1783ND	Knee High	15 - 20 mmHg	Nude	Diamond Pattern S-M-L-XL

LADIES' OPAQUE THERAPEUTIC SUPPORT HOSIERY PG. 10 - 13

MODEL	LENGTH	COMPRESSION	COLOR	SIZES
0371BG	Knee High	15 - 20 mmHg	Beige	Open Toe S-M-L-XL
0371BL	Knee High	15 - 20 mmHg	Black	Open Toe S-M-L-XL
0373BG	Knee High	15 - 20 mmHg	Beige	S-M-L-XL
0373BL	Knee High	15 - 20 mmHg	Black	S-M-L-XL
0374BG	Thigh High	15 - 20 mmHg	Beige	S-M-L-XL
0374BL	Thigh High	15 - 20 mmHg	Black	S-M-L-XL
0361BG	Knee High	20 - 30 mmHg	Beige	Open Toe S-M-L-XL
0361BL	Knee High	20 - 30 mmHg	Black	Open Toe S-M-L-XL
0363BG	Knee High	20 - 30 mmHg	Beige	S-M-L-XL
0363BL	Knee High	20 - 30 mmHg	Black	S-M-L-XL
0362BG	Knee High	20 - 30 mmHg	Beige	Open Toe S-M-L-XL
0362BL	Thigh High	20 - 30 mmHg	Black	Open Toe S-M-L-XL
0364BG	Thigh High	20 - 30 mmHg	Beige	S-M-L-XL
0364BL	Thigh High	20 - 30 mmHg	Black	S-M-L-XL
0374BG	Thigh High	15 - 20 mmHg	Beige	S-M-L-XL
0374BL	Thigh High	15 - 20 mmHg	Black	S-M-L-XL
0365BG	Pantyhose	20 - 30 mmHg	Beige	S-M-L-XL
0365BL	Pantyhose	20 - 30 mmHg	Black	S-M-L-XL
0375BG	Pantyhose	15 - 20 mmHg	Beige	S-M-L-XL
0375BL	Pantyhose	15 - 20 mmHg	Black	S-M-L-XL

TRUFORM TRUSHEER THERAPEUTIC SUPPORT HOSIERY PG. 6 - 9

MODEL	LENGTH	COMPRESSION	COLOR	SIZES
0253BG	Knee High	30 - 40 mmHg	Beige	S-M-L-XL
0253BL	Knee High	30 - 40 mmHg	Black	S-M-L-XL
0253TP	Knee High	30 - 40 mmHg	Taupe	S-M-L-XL
0254BG	Thigh High	30 - 40 mmHg	Beige	S-M-L-XL
0254BL	Thigh High	30 - 40 mmHg	Black	S-M-L-XL
0254TP	Thigh High	30 - 40 mmHg	Taupe	S-M-L-XL
0255BG	Pantyhose	30 - 40 mmHg	Beige	S-M-L-XL
0255BL	Pantyhose	30 - 40 mmHg	Black	S-M-L-XL
0255TP	Pantyhose	30 - 40 mmHg	Taupe	S-M-L-XL
0263BG	Knee High	20 - 30 mmHg	Beige	S-M-L-XL
0263BL	Knee High	20 - 30 mmHg	Black	S-M-L-XL
0263TP	Knee High	20 - 30 mmHg	Taupe	S-M-L-XL
0264BG	Thigh High	20 - 30 mmHg	Beige	S-M-L-XL
0264BL	Thigh High	20 - 30 mmHg	Black	S-M-L-XL
0264TP	Thigh High	20 - 30 mmHg	Taupe	S-M-L-XL
0265BG	Pantyhose	20 - 30 mmHg	Beige	S-M-L-XL
0265BL	Pantyhose	20 - 30 mmHg	Black	S-M-L-XL
0265TP	Pantyhose	20 - 30 mmHg	Taupe	S-M-L-XL
0267BG	Pantyhose	20 - 30 mmHg	Beige	Maternity S-M-L-XL
0267BL	Pantyhose	20 - 30 mmHg	Black	Maternity S-M-L-XL

TRUFORM LADIES' THERAPEUTIC SUPPORT SOCKS PG. 14 - 15 & 19

MODEL	LENGTH	COMPRESSION	COLOR	SIZES
1963BL	Knee High	10 - 20 mmHg	Black	Casual Style S-M-L-XL
1963NV	Knee High	10 - 20 mmHg	Navy	Casual Style S-M-L-XL
1963TN	Knee High	10 - 20 mmHg	Tan	Casual Style S-M-L-XL
1963WH	Knee High	10 - 20 mmHg	White	Casual Style S-M-L-XL
1973BL	Knee High	15 - 20 mmHg	Black	Rib Pattern S-M-L-XL
1973BN	Knee High	15 - 20 mmHg	Brown	Rib Pattern S-M-L-XL
1973NV	Knee High	15 - 20 mmHg	Navy	Rib Pattern S-M-L-XL
1973TN	Knee High	15 - 20 mmHg	Tan	Rib Pattern S-M-L-XL
1973WH	Knee High	15 - 20 mmHg	White	Rib Pattern S-M-L-XL
1975BL	Knee High	15 - 20 mmHg	Black	Cable Pattern S-M-L-XL
1975BN	Knee High	15 - 20 mmHg	Brown	Cable Pattern S-M-L-XL
1975NV	Knee High	15 - 20 mmHg	Navy	Cable Pattern S-M-L-XL
1975TN	Knee High	15 - 20 mmHg	Tan	Cable Pattern S-M-L-XL
1975WH	Knee High	15 - 20 mmHg	White	Cable Pattern S-M-L-XL
1976BL	Knee High	15 - 20 mmHg	Black	Diamond Pattern S-M-L-XL
1976BN	Knee High	15 - 20 mmHg	Brown	Diamond Pattern S-M-L-XL
1976NV	Knee High	15 - 20 mmHg	Navy	Diamond Pattern S-M-L-XL
1976TN	Knee High	15 - 20 mmHg	Tan	Diamond Pattern S-M-L-XL
1976WH	Knee High	15 - 20 mmHg	White	Diamond Pattern S-M-L-XL

TRUFORM MEN'S THERAPEUTIC SUPPORT SOCKS PG. 16 - 19

MODEL	LENGTH	COMPRESSION	COLOR	SIZES
1932BL	Crew Length	15 - 20 mmHg	Black	Athletic Style S-M-L-XL
1932WH	Crew Length	15 - 20 mmHg	White	Athletic Style S-M-L-XL
1932BN	Crew Length	15 - 20 mmHg	Brown	Athletic Style S-M-L-XL
1932NV	Crew Length	15 - 20 mmHg	Navy	Athletic Style S-M-L-XL
1932TN	Crew Length	15 - 20 mmHg	Tan	Athletic Style S-M-L-XL
1933BL	Knee High	15 - 20 mmHg	Black	Casual Style S-M-L-XL
1933BN	Knee High	15 - 20 mmHg	Brown	Casual Style S-M-L-XL
1933NV	Knee High	15 - 20 mmHg	Navy	Casual Style S-M-L-XL
1933TN	Knee High	15 - 20 mmHg	Tan	Casual Style S-M-L-XL
1933WH	Knee High	15 - 20 mmHg	White	Athletic Style S-M-L-XL
1934BL	Knee High	20 - 30 mmHg	Black	Athletic Style S-M-L-XL
1934WH	Knee High	20 - 30 mmHg	White	Athletic Style S-M-L-XL
1942BL	Knee High	8 - 15 mmHg	Black	Dress Style S-M-L-XL
1942BN	Knee High	8 - 15 mmHg	Brown	Dress Style S-M-L-XL
1942NV	Knee High	8 - 15 mmHg	Navy	Dress Style S-M-L-XL
1942TN	Knee High	8 - 15 mmHg	Tan	Dress Style S-M-L-XL
1943BL	Knee High	15 - 20 mmHg	Black	Dress Style S-M-L-XL
1943CH	Knee High	15 - 20 mmHg	Charcoal	Dress Style S-M-L-XL
1943BN	Knee High	15 - 20 mmHg	Brown	Dress Style S-M-L-XL
1943NV	Knee High	15 - 20 mmHg	Navy	Dress Style S-M-L-XL
1943TN	Knee High	15 - 20 mmHg	Tan	Dress Style S-M-L-XL
1943WH	Knee High	15 - 20 mmHg	White	Dress Style S-M-L-XL
1944BL	Knee High	20 - 30 mmHg	Black	Dress Style S-M-L-XL
1944BN	Knee High	20 - 30 mmHg	Brown	Dress Style S-M-L-XL
1944NV	Knee High	20 - 30 mmHg	Navy	Dress Style S-M-L-XL
1944TN	Knee High	20 - 30 mmHg	Tan	Dress Style S-M-L-XL
1945BL	Thigh High	20 - 30 mmHg	Black	Dress Style S-M-L-XL
1954BL	Knee High	30 - 40 mmHg	Black	Dress Style S-M-L-XL
1954BN	Knee High	30 - 40 mmHg	Brown	Dress Style S-M-L-XL
1954NV	Knee High	30 - 40 mmHg	Navy	Dress Style S-M-L-XL
1954TN	Knee High	30 - 40 mmHg	Tan	Dress Style S-M-L-XL

TRUFORM UNISEX TRAVEL SERIES SUPPORT SOCKS PG. 20 - 21

MODEL	LENGTH	COMPRESSION	COLOR	SIZES
1923BL	Knee High	15 - 20 mmHg	Black	S-M-L-XL
1923BN	Knee High	15 - 20 mmHg	Brown	S-M-L-XL
1923CH	Knee High	15 - 20 mmHg	Charcoal	S-M-L-XL
1923NV	Knee High	15 - 20 mmHg	Navy	S-M-L-XL

MEDICAL STOCKINGS 15 - 20 MMHG PG. 24 - 26				
MODEL	LENGTH	TOE STYLE	COLOR	SIZES
8875BG	Below Knee	Closed Toe	Beige	S-M-L-XL
8875BL	Below Knee	Closed Toe	Black	S-M-L-XL
0875BG	Below Knee	Open Toe	Beige	S-M-L-XL
0875BL	Below Knee	Open Toe	Black	S-M-L-XL

MEDICAL STOCKINGS 20 - 30 MMHG PG. 24 - 28				
MODEL	LENGTH	TOE STYLE	COLOR	SIZES
0864BG	Below Knee	Open Toe	Beige	Dot Top S-M-L-XL-2XL-3XL
0865BG	Below Knee	Open Toe	Beige	S-M-L-XL-2XL-3XL-4XL
0865BN	Below Knee	Open Toe	Brown	S-M-L-XL-2XL-3XL
0865BL	Below Knee	Open Toe	Black	S-M-L-XL-2XL-3XL-4XL
0865WH	Below Knee	Open Toe	White	S-M-L-XL-2XL-3XL
0865GR	Below Knee	Open Toe	Gray	S-M-L-XL-2XL-3XL
0865CH	Below Knee	Open Toe	Charcoal	S-M-L-XL-2XL-3XL
0865SGB	Below Knee	Open Toe	Beige	Short Length S-M-L-XL
0865SBL	Below Knee	Open Toe	Black	Short Length S-M-L-XL
0866BG	Thigh High	Open Toe	Beige	S-M-L-XL
0868BG	Thigh High	Open Toe	Beige	Dot Top S-M-L-XL-2XL
0868BL	Thigh High	Open Toe	Black	Dot Top S-M-L-XL-2XL
0868WH	Thigh High	Open Toe	White	Dot Top S-M-L-XL-2XL
8864BG	Below Knee	Closed Toe	Beige	Dot Top S-M-L-XL-2XL-3XL
8865BG	Below Knee	Closed Toe	Beige	S-M-L-XL-2XL-3XL-4XL
8865SGB	Below Knee	Closed Toe	Beige	Short Length S-M-L-XL
8865SBL	Below Knee	Closed Toe	Black	Short Length S-M-L-XL
8865BL	Below Knee	Closed Toe	Black	S-M-L-XL-2XL-3XL-4XL
8865WH	Below Knee	Closed Toe	White	S-M-L-XL-2XL-3XL
8865BN	Below Knee	Closed Toe	Brown	S-M-L-XL-2XL-3XL
8865CH	Below Knee	Closed Toe	Charcoal	S-M-L-XL-2XL-3XL
8865GR	Below Knee	Closed Toe	Gray	S-M-L-XL-2XL-3XL
8866BG	Thigh High	Closed Toe	Beige	S-M-L-XL
8867BG	Thigh High	Closed Toe	Beige	Lace Top S-M-L-XL
8868BG	Thigh High	Closed Toe	Beige	Dot Top S-M-L-XL-2XL
8868BL	Thigh High	Closed Toe	Black	Dot Top S-M-L-XL-2XL
8868WH	Thigh High	Closed Toe	White	Dot Top S-M-L-XL-2XL

MEDICAL STOCKINGS 30 - 40 MMHG PG. 24 - 28				
MODEL	LENGTH	TOE STYLE	COLOR	SIZES
0844BG	Below Knee	Open Toe	Beige	Dot Top S-M-L-XL-2XL-3XL
0845BG	Below Knee	Open Toe	Beige	S-M-L-XL-2XL-3XL
0845BL	Below Knee	Open Toe	Black	S-M-L-XL-2XL-3XL
0845WH	Below Knee	Open Toe	White	S-M-L-XL-2XL-3XL
0845SGB	Below Knee	Open Toe	Beige	Short Length S-M-L-XL
0846BG	Thigh High	Open Toe	Beige	S-M-L-XL
0848BG	Thigh High	Open Toe	Beige	Dot Top S-M-L-XL
8844BG	Below Knee	Closed Toe	Beige	Dot Top S-M-L-XL-2XL-3XL
8845BG	Below Knee	Closed Toe	Beige	S-M-L-XL-2XL-3XL
8845SGB	Below Knee	Closed Toe	Beige	Short Length S-M-L-XL
8845BL	Below Knee	Closed Toe	Black	S-M-L-XL-2XL-3XL
8845WH	Below Knee	Closed Toe	White	S-M-L-XL-2XL-3XL
8846BG	Thigh High	Closed Toe	Beige	S-M-L-XL
8848BG	Thigh High	Closed Toe	Beige	Dot Top S-M-L-XL
8848BL	Thigh High	Closed Toe	Black	Dot Top S-M-L-XL

MEDICAL PANTYHOSE 20 - 30 MMHG PG. 29				
MODEL	LENGTH	TOE STYLE	COLOR	SIZES
1756BG	Pantyhose	Closed Toe	Beige	P-M-T-XT
1756BL	Pantyhose	Closed Toe	Black	P-M-T-XT
1757BG	Pantyhose	Closed Toe	Beige	Maternity P-M-T-XT
1757BL	Pantyhose	Closed Toe	Black	Maternity P-M-T-XT
1758BG	Pantyhose	Closed Toe	Beige	Full Figure P-M-T-XT
1758BL	Pantyhose	Closed Toe	Black	Full Figure P-M-T-XT

ANTI-EMBOLISM STOCKINGS 20 - 30 MMHG PG. 30 & 31				
MODEL	LENGTH	TOE STYLE	COLOR	SIZES
0808BG	Below Knee	Open Toe	Beige	S-M-L-XL-2XL-3XL
0808WH	Below Knee	Open Toe	White	S-M-L-XL-2XL-3XL
0808BGS	Below Knee	Open Toe	Beige	Short Length S-M-L-XL
0810BG	Thigh Length	Open Toe	Beige	S-M-L-XL
0810BL	Thigh Length	Open Toe	Black	S-M-L-XL
0810WH	Thigh Length	Open Toe	White	S-M-L-XL
8808BG	Below Knee	Closed Toe	Beige	S-M-L-XL-2XL-3XL-4XL
8808BL	Below Knee	Closed Toe	Black	S-M-L-XL-2XL-3XL-4XL
8808WH	Below Knee	Closed Toe	White	S-M-L-XL-2XL-3XL-4XL
8808BGS	Below Knee	Closed Toe	Beige	Short Length S-M-L-XL
8810BG	Thigh Length	Closed Toe	Beige	S-M-L-XL
8810BL	Thigh Length	Closed Toe	Black	S-M-L-XL
8810WH	Thigh Length	Closed Toe	White	S-M-L-XL

NEW! PREMIUM MICROFIBER MEDICAL SOCKS 20 - 30 MMHG PG. 22 & 23				
MODEL	LENGTH	TOE STYLE	COLOR	SIZES
9965BG	Below Knee	Closed Toe	Beige	S-M-L-XL-2XL-3XL
9965BL	Below Knee	Closed Toe	Black	S-M-L-XL-2XL-3XL
0965BG	Below Knee	Open Toe	Beige	S-M-L-XL-2XL-3XL
0965BL	Below Knee	Open Toe	Black	S-M-L-XL-2XL-3XL

HEALTH & DIABETIC TRUSOFT & COMFORT CARE SOCKS PG. 32 - 34				
MODEL	LENGTH	COMPRESSION	COLOR	SIZES
1911WH	Mini-Crew	N/A	White	Trusoft XS-S-M-L-XL
1912BL	Crew Length	8 - 15 mmHg	Black	Trusoft XS-S-M-L-XL
1912WH	Crew Length	8 - 15 mmHg	White	Trusoft XS-S-M-L-XL
1913BL	Knee High	8 - 15 mmHg	Black	Trusoft XS-S-M-L-XL
1913WH	Knee High	8 - 15 mmHg	White	Trusoft XS-S-M-L-XL
1918BL	Crew Length	N/A	Black	Loose Fit Top S-M-L-XL
1918WH	Crew Length	N/A	White	Loose Fit Top S-M-L-XL
1918TN	Crew Length	N/A	Tan	Loose Fit Top S-M-L-XL
1918NV	Crew Length	N/A	Navy	Loose Fit Top S-M-L-XL

READY-WEAR UPPER EXTREMITY PG. 35				
MODEL	LENGTH	COMPRESSION	COLOR	SIZES
3315BG	Arm Sleeve	15 - 20 mmHg	Beige	S-M-L
3315BL	Arm Sleeve	15 - 20 mmHg	Black	S-M-L
3315BN	Arm Sleeve	15 - 20 mmHg	Brown	S-M-L
3316BG	Arm Sleeve	15 - 20 mmHg	Beige	Dot Top S-M-L
3325BG	Arm Sleeve	20 - 30 mmHg	Beige	S-M-L
3325BL	Arm Sleeve	20 - 30 mmHg	Black	S-M-L
3325BN	Arm Sleeve	20 - 30 mmHg	Brown	S-M-L
3326BG	Arm Sleeve	20 - 30 mmHg	Beige	Dot Top S-M-L
3326BL	Arm Sleeve	20 - 30 mmHg	Black	Dot Top S-M-L
3326BN	Arm Sleeve	20 - 30 mmHg	Brown	Dot Top S-M-L
3312BG	Guantlet	15 - 20 mmHg	Beige	S-M-L
3322BG	Guantlet	20 - 30 mmHg	Beige	S-M-L
3322BL	Guantlet	20 - 30 mmHg	Black	S-M-L
3322BN	Guantlet	20 - 30 mmHg	Brown	S-M-L

ACCESSORIES PG. 36		
MODEL	DESCRIPTION	SIZES
0750	Slip Sock	N/A
0752	Stocking Donner	R-T-L
0757	Donning Gloves	S-M-L-XL
0758	It Stays! Roll-On Body Adhesive	N/A
0760	3" Garter Belt	XS-S-M-M/L-L-XL
0762	2" Lightweight Garter Belt	S-M-L-XL
0765	6" Garter Belt	XS-S-M-M/L-L-XL



COMPRESSION
THERAPY COLLECTION



TRUFORM is a Division of SAI Therapeutic Brands
ph: 800-888-0458 • fax: 800-309-9055 • truform.com

3960 Rosslyn Dr, Cincinnati, OH
45209, United States

In Canada, a Division of Airway Surgical Appliances Ltd.
ph: 800-267-3476 • fax: 800-267-7805 • truform.com

189 Colonnade Rd S, Ottawa, ON
K2E 7J4, Canada



TRUFORM[®]
LEG HEALTH